





HEALEON
MEDICAL

Leadership in Regenerative Solutions





CENTRIfyte™ 1000

This all-in-one device simplifies the isolation of cells from all tissue sources using a semi-automatic methodology, which eliminates the need for additional personnel and allows for simple, consistent reproducible results. The **CENTRIfyte™** has been proven to provide consistent higher yields compared to that of competitors from the same sample volume. The **CENTRIfyte™**, when coupled with the expertise of Healeon Medical, enables the clinician to maximize cell counts and viability from the smallest of harvested tissue samples.

Tissue Sources:

- Bone Marrow Concentrate
- Stromal Vascular Fraction
- Platelet Rich Plasma
- Adipose Tissue Concentrate

The simplified user interface allows for easy navigation through associated steps. Each program has been preset to meet the specific protocol standards. Simply choose the specimen protocol needed and a simple push of the button processes all autologous needs.

Proprietary specimen kits allow ease of transfer through a complete closed system, which minimizes cross contamination as well as reduction of cell viability due to environmental influences.



HEALEON MEDICAL PRP

Healeon Medical's PRP System is a proprietary platform based on blood separation according to density, using a unique separation gel. The system allows for the preparation of pure PRP with optimal platelet concentrations and cellular content. Virtually, all Red Blood Cells (RBC) are eliminated from PRP and the concentration of pro-inflammatory Granulocytes (GRs) is dramatically reduced.

PRP is prepared by taking a small sample of the patient's own blood, then separating platelets from Platelet-Poor Plasma (PPP), RBC and GRs via centrifugation. PRP is then collected and can be applied the treated site to promote soft and hard tissue regeneration. The whole preparation process takes up to 10 minutes.

EFFECTIVENESS:

- Optimum platelet concentration of 2.5-5x baseline - only one centrifugation cycle
- Rapid and simple process - only one centrifugation
- Adjustment to a specific clinical application by controlling the final PRP volume
- Compatible with OR and/or doctor's office setting

HEALEON MEDICAL ACCESSORIES



AUTOMATED CELL COUNTER:

Incredible performance in a small package.

- More accurate and reliable than any automated cell counter
- Results in just 8 seconds
- Ultra-compact, ultra-easy to use (no manual counting, no focusing of any kind)
- Complete, high resolution cell size/count histogram
- Post-processing/analysis of count distributions with user-adjustable regions/gates
- Complete data analysis & data portability (connectivity to PC and MoxiChart software package)
- No maintenance required



HEALEON PROCEDURE PACK:

Customized surgical procedure pack for the minor in-office surgery. This all-in-one pack eliminates the need to source multiple sterile materials - it is all included in an inexpensive set. *Simply open, prep, drape and treat.*

SPECIMEN HARVESTING and REINJECTION ACCESSORIES:

Healeon Medical can provide harvesting cannulas and reinjection needles made to your specifications.

"As a Board Certified Plastic Surgeon, it has been very exciting to see the advent of the minimally invasive, regenerative therapy modalities, whether it be the PRP or fat-derived SVF's using the CENTRicyte. The applications are very diverse and it has been a fantastic addition to our surgical practice."

Dr. Semira Bayati, Newport Beach, CA



**HEALEON
MEDICAL**

Website: www.Healeonmedical.com

"I have been a client of Healeon Medical since 2012. They have been instrumental in keeping my Regenerative Medicine/Cosmetic practice a leader in the field. Healeon always keeps us informed about emerging technology in the field. They are quick to respond whenever there are questions or concerns. They take a personal interest in the success of their clients.

I recently purchased a CENTRicyte from Healeon, which puts me at the forefront of regenerative medicine technology. They have been outstanding with customer support and training. The transition to the CENTRicyte was effortless due to the efforts of the Healeon team. I would recommend Healeon and the CENTRicyte to any of my colleagues".

Dr. John P. Laura, Newport Beach, CA

1111 Rancho Conejo Blvd., Suite 204
Newbury Park, CA 91320
Phone: 805-499-7482

Email: sales@healeonmedical.com